

**medima**

Safe and precise infusion

*Medima Line Sets*



Administration of drugs, liquids and blood  
Parenteral and enteral feeding  
Oncological treatment

# Medima Line administration sets for use with volumetric infusion pumps

## Standard sets

MEDIMALINES  
PVC DEHP-free, 270cm  
P.P2 REF: 0800100  
P100, P200, P300 REF: 08001100

MEDIMALINE L  
PVC DEHP-free, Light protection  
270cm, P.P2 REF: 08003000  
P100, P200, P300 REF: 08003100

ST110  
PVC DEHP-free, 285cm  
REF: 08110000

ST110L\*  
PVC DEHP-free,  
Light protection, 285cm,  
REF: 08110100

ST110P\*\*  
PVC-free, 285cm  
REF: 08110200

ST110PL\*\*  
PVC-free, Light protection, 285cm  
REF: 08110300


ST111  
PVC DEHP-free, 285cm  
REF: 08111000

ST111L\*  
PVC DEHP-free,  
Light protection, 285cm  
REF: 08111100

ST112  
PVC DEHP-free, 285cm  
REF: 08111200

ST112P\*\*  
PVC-free, needle-free, 285cm  
REF: 08111220

ST112L  
PVC DEHP-free, Light protection  
needle-free, 285cm  
REF: 081112100



**Blood set**

MEDIMA LINE B  
PVC DEHP-free, 270cm  
P100, P22 REF: 08002000  
P100, P200, P300 REF: 08002100

**Oncological sets**

ON25 PVC DEHP-free, needle-free, 304cm  
REF: 083325000

ON25P\*\*  
PVC-free, needle-free, 304cm  
REF: 083325200

ON25L\*  
PVC DEHP-free, needle-free, 304cm  
REF: 083325100

ON25PL\*\*  
PVC-free, Light protection  
needle-free, 304cm  
REF: 083325300

**Nutrition set**

NE52\*  
PVC DEHP-free,  
needle-free (ENFit), 285cm  
REF: 085520000


**Extension lines**

EX21 PVC DEHP-free, needle-free, 40cm  
REF: 08921000

EX21L\*  
PVC DEHP-free, Light protection  
needle-free, 40cm  
REF: 08921100

EX21P\*\*  
PVC-free, needle-free, 40cm  
REF: 08921200

EX21PL\*\*  
PVC-free, Light protection  
needle-free, 40cm  
REF: 08921300



Medima Line administration sets  
are equipped with:

- Free flow protection clamp (FFPC)
- Roller clamp with spike handler
- Luer Lock cap with hydrophobic filter  
(automatic prime stop function)\*\*\*
- Rotating Luer Lock\*\*\*

**medima**

**medima**

Medima Sp. z o.o.  
Al. Jerozolimskie 200, 02-486 Warsaw, Poland  
Phone: +48 (22) 313 22 66 Fax: +48 (22) 313 22 69  
email: medima@medima.pl  
[www.medima.pl](http://www.medima.pl)

## Volumetric Infusion Pumps

Medima Line administration sets are compatible with all models of volumetric infusion pumps by Medima



### P100

- Continuous infusion in ml/h
- Programmable flow rate based upon volume and time
- Flow rate range 0.1 – 1200 ml/h

### P200

- P100 capabilities
- Infusion in all popular units
- Infusion modes: continuous, intermittent, profile, TPN

### P300

- P200 capabilities
- Drug library (DERS)



### P

- Infusion in ml/h
- Flow rates 0.1-1200 ml/h

### P2

- P capabilities
- Infusion in all popular units
- Infusion modes: continuous, intermittent, profile