


Safe and precise infusion

*MEDIMA Line IV Sets*



**medima**



### MEDIMA Line S

- Standard IV set for infusion of drugs and liquids
- Drop chamber, 20 drops/ml, filter 15 $\mu$ m
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free



### MEDIMA Line B

- IV set for infusion of blood and blood similar liquids
- Drop chamber, 20 drops/ml, filter 200 $\mu$ m
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free



### MEDIMA Line L

- IV set for infusion of light sensitive drugs, amber
- Drop chamber, 20 drops/ml, filter 15 $\mu$ m
- Free flow protection clamp
- Roller clamp
- Luer lock
- Infusion port
- Length: 270 cm
- DEHP free
- Latex free

MEDIMA Line IV sets can be used for administration of cytostatics

Medima Sp. z o.o.

Al. Jerozolimskie 200, 02-486 Warsaw, Poland  
tel. +48 (22) 313 22 66; fax +48 (22) 313 22 69  
e-mail: [medima@medima.pl](mailto:medima@medima.pl); <http://www.medima.pl>



**INNOVATIVE ECONOMY**  
NATIONAL COHESION STRATEGY



EUROPEAN UNION  
EUROPEAN REGIONAL  
DEVELOPMENT FUND



Project co-financed by the European Union from European Regional Development Fund